

PMU

Permenent Make-Up Devices

Beginners

GMD-S: Standard type for beginner

- Handy & portable to use
- 3 steps of speed
- 4500 ~ 6000 RPM



Professional

GMD-H: Handpiece type with dual functions

- Avaliable for Micropigmentation and Microneedling
- 5 steps of speed
- 4000 ~ 6400 RPM



GMD-T: Triple functions & 2.5mm needle length

- Micropigmentation
- Microneedling
- Tatto



GMD-R: Smart memory function & User friendly interface


- Full touch screen operation
- Avaliable upto 2.5mm needle depth & 13 steps speed
- 3 mode : Easy mode, Professional mode and User mode
- Medical use





PMU

3 ADVANTAGES OF DMD DEVICE




01 Fastening Method using Magnet

- Needle and Device is connected by Magnetic
- Unshaking during operation
- Strong motor power can be delivered to the needles



02 PID(Proportional integral derivative) control

- PID controls maintains the same power and speed under all conditions that the pigments penetrates the skin evenly



03 Adjustable needle depth

- Adjustable needle depth by 0.1mm upto 2.0mm
- Tick marks in 0.2mm

VARIOUS TYPES OF TIPS

	WTNR-01 WTNR-01N		NN WTN	N WTNR-01EN		IR-03	WTNR-03EN		WTNR-05	
ROUND	•	•		•		8	*		*	
	0.4mm/27G	0.2mm/330	G 0.151	0.15mm/35G		m/31G	0.15mm/35G		0.25mm/31G	
	Eyeliner, Eyebrow	Feathering Brow, Fine hair/ Eyelin	e tre	No-pain treatment Eyebrow		liner, eliner	No-pain treatment Lip		Lip, Eyeliner	
FLAT / MAGNUM	WTNR-04F WTNF		R-04S	04S WTNR		F WTNR-06A		WTNR-09M		
	••••		•••		•••			****		
	0.3mm/30G	0.3mr	m/30G	0.25m	m/31G	0.2	0.25mm/31G		0.25mm/31G	
	Shadowing Brows Hairline, Lip-coloring	EMBO	Hai		ng Brows ie, Lip, Tattoo	E	EmboBrows, Hairline, Lip		Lip, Areola, Tattoo	
SPECIAL	WTNR-07		WTNR-22A			WTNR-09(2.0)		MASSAGE		
						• • •				
	0.25mm/31G		0.25mm/31G			0.2mm/33G		Silicon		
	EmboBrows, Hairline, Lip		Tattoo, Lip and Areola			Microneedling		Skin Massage		