



Woorhi Mechatronics Co., Ltd. is a global manufacturer of microneedling, micro-pigmentation products. Established in 2008, Woorhi Mechatronics is the pathfinder in the medical and aesthetic market.

Woorhi Mechatronics have state of the art facilities with ISO 13485, GMP Certified.

Currently, we are exporting our products to over 100 companies in 30 countries such as the United States, United Kingdom, Austrailia and providing customized follow-up service.



Clean Room



Assemble Line



CERTIFICATE AND AWARDS

- ISO13485, 14001 and 9001
- MFDS, CE medical, CPNP and Dermatest
- Prime Minister's Commendation



CERTIFICATE OF BRAND K

ISO 13485



ISO 14001



Inno - Biz Certificate



ISO 9001



MICRONEEDLING

01 RAFFINE

02 CODE-X

03BELLA

04 MINIPIN

05 NOBLE

PERMANENT MAKE UP

01 DMD-S

02 DMD-H

03 DMD-T

04 DMD-R

COSMETICS AND ACCESSORIES

01 AMPOULE

O2 SERUM

03 CREAM

04 FACIAL SHEET MASK

05 SCALP CARE

06 MICRONEEDLING ACCESSORIES







1 MICRONEEDLING For PROFESSIONAL



O1 RAFFINE

Medical Grade Microneedling Pen

- One of the few CE marked microneedling devices in the market(CE MDD 93/42/EEC)
- Plug-in system for full-power treatments and consistent output
- High quality motor that guarantees smooth skin needling even with the deepest needle depth
- 0 2.0mm adjustable needle depth

C€ 0068



02 **CODE-X**

New & Innovative Wireless Microneedling Pen

- High capacity lithium ion battery(600mAh, 2.5H run time)
- Battery Management System that maintain same power regardless of remaining battery level
- Ergonomic and lightweight design
- High quality motor that guarantees smooth skin needling even with the deepest needle depth
- 0 2.0mm adjustable needle depth



03 **BELLA**

Professional Microneedling Device for Aestheticians

- Ergonomic design and chromium plating
- High quality motor that guarantees smooth skin needling even with the deepest needle depth
- 0 2.0mm adjustable needle depth

AUTOMATIC MICRONEEDLING TREATMENT



Microneedling Treatment simulates skin cell and induce collagen. Thousand of tiny channels are made in a short time.

Through powerful motors, no damage on the skin and just penetrate skin layer from epidermis to deeper with a cosmetic and medical solutions.

Purpose of Microneedling Treatment

Skin Care: Skin Rejuvenating, Moisturizing, Whitening Trouble Care : Pigment. Acne Scar, Pimples

Body & Scalp Care: Cellulite Reduction, Hair-Growth

Before and After









BEFORE

AFTER

AFTER









BEFORE

AFTER

1 AT-HOME MICRONEEDLING DEVICES



Innovative All-In-One microneedling device

- All-in-one design needle cartridges that maximizes user convenience
- 12pin microneedles(fixed at 0.25mm)
- Cosmetic delivery system approved by clinical results.
- Patented new operating system using solenoid electromagnet




■ Technical data of Microneedling device

	RAFFINE	CODE-X	BELLA	MINIPIN	NOBLE	
Features	Medical use	Aesthe	tic, Spa	Home use		
Needle Adjustment	Ad	justable up to 2.0r	0.25mm	0.2mm		
Certification	KC, CE MDD	KC, CE MFDS, FDA(2023) CE MDR(2024)		KC, CE		
Speed(RMP) (Rotating per Minute)	4500 – 6500	3500-5500	5000-6500	3500	4500-6500	
Weight	80g	65g	127g	60g	76g	
Wire/Wireless	Wire	Wireless	Wire	Wireless	Wire	
Hand-piece Size(cm)		15.5	6	15.5		
Power Input		5VDC	3.7VDC	5VDC		

Guideline for Microneedling

	Treatment Area	Recommendation	Duration	Length	
	Moisturizing	3 times or more	Every 1 weeks	0.2-0.6mm	
Basic Care	Whitening	3 times or more	Every 1 weeks	0.2-0.6mm	
	Wrinkles	5 times or more	Every 1 weeks	0.2-1.2mm	
	Dark Circles	5 times or more	Every 1 weeks	0.2-0.6mm	
	Skin Rejuvenating	5 times or more	Every 2 weeks	0.2-1.5mm	
Intensive care	Alopecia	10 times or more	Every 1 weeks	0.2-1.0mm	
	Hyperpigmentation	10 times or more	Every 1 weeks	1.5-2.0mm	
	Stretch Mark	10 times or more	Every 2 weeks	1.5-2.0mm	
	Cellulite Reduction	10 times or more	Every 2 weeks	1.5-2.0mm	
	Scars	10 times or more	Every 2 weeks	1.5-2.0mm	



2 PMU



01 **DMD-S**

Standard type for beginner

- Handy & portable to use
- 3steps of speed
- on/off button type to easy control



02 **DMD-H**

Handpiece type with dual functions

- 5 steps of speed
- On/off button type to easy control
- Available for Micropigmentation, Microneedling



03 **DMD-T**

Triple functions & 2.5mm needle length

- Micropigmentation
- Microneedling
- Tatto

04 **DMD-R**

Smart memory function & friendly user interface

- Full touch screen operation
- Available upto 2.5mm needle depth & 13 steps speed
- Medical Use









■ What is DMD (Digital Micropigmentation Device)?

Digital Micropigmentation Device makes the hole and leave the pigments(ink) into the skin layer. The pigments(color) will be stayed at least 6 to 36months.

The period is flexible following the skin condition of customers.

As a digital type of the device, it can control the speed, power, and needle depth during the treatment. DMD have dual function of Microneedling and Permanent make up device.

Special Technology of Device

DMD has the special technology called PID control.

PID(Proportional integral derivative control) can adjust and stay the movement of motor. It is the ideally products for the permanent make up, because that if the motor cannot control by program, the speed and power of needles is uneven and then the color is not stable too. Also we developed the different type of setting part of the needles, and it could be made the treatment without any problem as like pop-out or dragging on the skin.

Before and After

Eyebrow









AFTER

Eyeline





BEFORE

AFTER

2 PMU

■ Technical data of Microneedling device

	DMD-S	DMD-H DMD-T		DMD-R	
Features	Beginner	Beginner Professional	Professional Tatto treatment		
Needle Adjustment	Adjustable (up to 2.0mm	Adjustable up to 2.5mm		
Certification	CE, N	MFDS	CE	CE, MFDS	
Speed(RMP) (Rotating per Minute)	4500 - 6000	4000-6400	4000-6500	3000 – 6500	
Weight(Pen)	66g	64g	72g	82g	
Hand-piece Size(cm)		13.5		14	
Power Input	6V	DC	24VDC	6VDC	
Specialized	Handy, 3steps speed	Smart Program, Handy	2 pcs Handpiece Smart Program	Smart Program, Memory of Setting	

■ Guideline for Permanent make up

TYPE OF DISPOSABLE TIPS

ROUND Eyeliner, Eyebrow	WTNR-01NN		WTNR-01EN		WTNR-03		WTNR-03EN		WTNR-05	
	•			•		*				*
	-	Feathering Brow, Fine hair/ Eyeliner		No-pain treatment Eyebrow		Lip liner, Eyeliner		No-pain treatment Lip		Lip, Eyeliner
FLAT / MAGNUM Shadowing Brows, Hairlin Lip-coloring	WTNR-04F W		TNR-04	NR-04S WTNR-		06F WTNR-06A		NR-06A	WTNR-09M	
	••••		••••		••					
	EmboBrows Hairline, Lip		Hairline I II		Lip,	ip, Embobrows,		Lip, Areola, Tattoo		
	WTNR-0		7 WT		TNR-22A		WTNR-09(2.0)		MASSAGE	
SPECIAL										
EmboBr Hairline					too, d Areola Mic		Microneedling		Skin Massage	

Protocols for Permanent make up



Consultation with customers and design



STEP 2
Cleansing
Process



Anesthetic Cream



STEP 4
Treatment



Applying the recovery cream



Touch up treatment (1 month later)


PIGMENTS LINE UP

Organic and Long Lasting Color

Got safety test for all colors and ingredients.

11 basic colors and other special colors are available for customer's request.

Packed in 5ml capacity, it can be used hygienically without worrying about contamination after opening.



Don't need to mix the color and with only one color is showed well-colored what you want. From beginner to professionals, AWY Pigments provide the satisfied result!

AWY PIGMENT COLOR CHART



3 COSMETICS





MICRONEEDLING AMPOULE

[RealTox Anti-Wrinkle Ampoule] skin booster, wrinkle care, skin pore

The main ingredient is 'Real Tox', a patented ingredient extracted from plants, which has the effect of improving skin elasticity and wrinkles. Niacinamide, Protease and Bifida Ferment Lysate is included for anti-oxidant.

*Clinical report for skin-pore care

USE FOR

- Systematic management with Microneedling treatment
- Moisturizing
- Whitening
- Anti-wrinkle and Anti-aging
- Skin-elasticity
- Rejuvenating and Skin-texture care



02

MICRONEEDLING AMPOULE

[Super Hyalu Energy]

moisturizing and rejuvenating

Glyceryl Glucoside opens the moisture path in the skin. Pink Flower Biome fermented with a special method has anti-inflammatory, anti-oxidant and anti-wrinkle properties.

Lactobacillus, exists more in healthy skin in their 20s, improves skin cell regeneration, wrinkle improvement.

*Certificate: CPNP, Dermatest

USE FOR

- Systematic management with Microneedling treatment
- Moisturizing
- Rejuvenating
- Anti-wrinkle and Anti-aging
- Skin-elasticity
- Skin-texture care

03

PRE/POST CARE AND DAILY SERUM

[Collagen Lifting Facial Serum]

skin-lifting, rejuvenating and anti-aging

The **collagen water**, which is 80.5% of the total contents, replaced purified water and collagen invigorates dry and thin skin and protect the skin layer. **13 kinds of peptides** are added to maintain young and strengthen skin regenerative power. In addition, **9 plant stem cells** prevent skin aging and make the skin surface healthy with Centella Asiatica Extract that increases absorption rate by exosomes.

*Certificate: CPNP. Dermatest

USE FOR

- Elasticity and rejuvenating
- Moisturizing
- Anti-wrinkle and Anti-aging





04

PRE/POST CARE AND DAILY SERUM

[Vitamin Tree Brightening Serum] brightening, dark spot care and wrinkle care

Contains 78.9% of **vitamin tree extract** with strong antioxidant function, replacing purified water to help skin regeneration, anti-aging, and wrinkle relief. In addition, **triple brightening liposome** contains niacinamide and arbutin and cabbage fermented water containing a large amount of Lactobacillus lactic acid bacteria, vitamins, and various nutrients.

*Certificate: CPNP, Dermatest

USE FOR

- Brightening and Skin-tone care
- Pigmentation removal
- Moisturizing
- Anti-wrinkle and Anti-aging
- Skin-elasticity
- Skin-texture care

3 COSMETICS





05

POST CARE AND DAILY CREAM

[CICA Water Derma Cream]

Skin calming, Moisturizing

Cica Biome, which is effective in suppressing skin inflammation and recovering skin wounds.

Ceramide complex reconstructs damaged skin barriers. **Panthenol** maintains skin moisture and beneficial bacteria such as Bifida fermentation filtration.

*Certificate: CPNP, Dermatest

USE FOR

- Soothe the skin
- Dry skin care
- Moisturizing
- Wounds care

06

POST CARE AND DAILY MASK SHEETS

[CELLUVEN FACIAL MASK]

Skin cooling, Soothing effects

A mask sheet coated with **bio-cellulose** that improves adhension to skin and delivery system of active ingredients. Compared to Non-woven mask sheet, Celluven minimizes skin irritation and prevention of moisture evaporation.

*Certificate: OTC

USE FOR

- Soothe the skin
- Moisturizing
- Sensitive skin care with safety ingredients





ALWAYS



07



SCALP CARE SOLUTION

[Cytokines-8] POWDER GROWTH FACTOR

The main ingredient is 8 various cytokines for helping the hair growth and strengthening the hair roots. (KGF, FGF9, IGF-1, VEGF, Bfgf, Noggin, SOD and ATP) It is effective for male pattern hair loss due to aging, stress, and heredity, and the deeper it is absorbed into the scalp, the better the hair loss management effect.

ALWAYS





A tonic and spray that soothes sensitive scalp and relieves hair loss symptoms. Containing 22 natural ingredients, it improves the scalp to a healthy condition and provides a refreshing feeling the moment you apply it. Contains L-menthol, salicylic acid, and dexpentanol that are effective for hair loss management.

Before and After

MALE




BEFORE

BEFORE




AFTER





AFTER

FEMALE





AFTER





BEFORE

AFTER

14

4 AESTHETIC, SPA CARE



$\bigcirc 1$

SCALP DAILY CARE SYSTEM

Cooling and Soothing for Scalp care

This product helps soothing the scalp area. Heat-irritated scalp is known to cause hair loss. It is a product that immediately lowers the heat of the irritated scalp and keeps it cool.

Specification data sheet				
Recommended user	Personal / Clinical use			
Use time	10mins			
Electric type	Battery and requested charging			
Weight	500g			





02

PROTECTIVE SLEEVE

Safety sleeve/sheath for microneedling pen

This protective sleeve is designed for microneedling pens and PMU machines. It covers the entire pen and disposable tip in order to form a protective barrier against cross-contamination (no contact with blood and other fluids during the treatment)

Specification data sheet					
Package	200pcs in one box				
	Single use only				



03

SKIN COOLING STICK

Cooling and Soothing for skin and body care

AWY Skin Cooling Stick is used to soothe the skin.

After MTS or skin care, use COOLING STICK on your face for about 10 minutes to relieve swelling or heat. If placed in the freezer for 20 minutes before use, it will remain cold for up to 15 minutes after removal.

Specification data sheet				
Storage	freezer about 20mins before use			
Material	Aluminum(no heavy metals)			
Size/Weight	55mm x 95mm / 75g			



04

PORE CLEANSING FACIAL SPONGE

AWY

Cooling Stick

Daily facial cleanser with organic konjac sponge

Made from 100% natural konjac, this product is gentle on the skin and can be used daily. It does not contain any chemical additives, so it can be used even for sensitive skin. It can be used in a variety of ways depending on the purpose of use and skin condition as it is made of soft, rough and both sides. (Normal cotton/Walnut side)

Specification data sheet				
Wet thoroughly with lukewarm water before use.				
Material	100% natural organic konjac			
Size/Weight	8cm x 10cm			
Recommend	Replace every 2months			
	·			









Woorhi Mechatronics Co., Ltd.

5th Floor, Woorhi Mecha building, 58, Sanbon-ro, 48beon-gil, Gunpo-si, Gyeonggi-do, 15846, Republic of Korea

♦ +82-31-457-5032 **№** saleswrmc@gmail.com



International Sales Team

C +82-10-5446-5032

🖄 saleswrmc@gmail.com



Domestic Sales Team

C +82-10-8774-5032

🖄 awy1204@naver.com



Global Website



Domestic Brand Website



Company Introduction Video



Instagram

